

Visor For Medical, Odontological and Veterinary Mask ref. 22201



Replacement for medical, dental, and veterinary skullcaps

Materials

Made of PETG, it allows to be light, flexible and resistant to splatters of fluids generated during its use. Neutral color (Transparent), optically rectified, with protective film on both sides of the visor, guarantees total visibility without distortion during operation.

Parameter

Traction resistance to deformation MPa 53. Traction resistance to break point MPa 26. Elongation to break point % 200. Elasticity module in traction MPa 2,200. Resistance to flexing point MPa 79. Charpy impact resistance with notch kJ/m2 10. Rockwell hardness, scale M/R M/R 115

Maintenance

Must be washed with cold water and detergent-free soaps, after each use Air dry. Do not clean the visor dry or use petroleum derivatives such as gasoline or thinner for cleaning.

Normative compliance

ANSI/ISEA Z87.1-2015- Face Protection Devices - Full Requirement ANSI/ISEA Z87.1-2015- Face Protection Devices -Droplets and Splashes OSHA Regulation 29 CFR 1910 Subpart I 29 CFR 1910.133.