

## Face Mask ref. 22501



The face mask should be put on before wear the gloves, and should be removed after the gloves are removed. They should be changed between patients or when visibly stained. The shelf life of each mask is estimated at maximum use of 8 hours.

- **Materials**

Non-woven fabric 100% polypropylene  
Outer layer 25g/m<sup>2</sup> (±2%) Inner layer 25g/m<sup>2</sup> (±2%) Fabric filtration 25g/m<sup>2</sup> (±2%) Size 17.7 cm x 9.5 x 9 x 17.5 cm  
Number of strips 2 elastic units Length after elastic 17.5 ± 2cms. Material name Spun bond 100% POLYPROPYLENE (PP) Minimum nasal bridge length 8 cm, maximum 11 cm. Point breaking force between not less than 10N Anti-moisture. Mask particle filtration efficiency > 98%

- **Maintenance**

Device for single use, do not store in humid environments, if the packaging is open do not use. Discard the product after 8 hours of use.

- **Normative Compliance**

Health Record: 20101747

Invima Sanitary Registration: INVIMA 2015DM-0013949

Filtration level: ≥ 98%

Dimensions: 17,5 x 9,5 Cm ± 2mm